



Patient Safety

- Reduces the risk of CO₂ rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather
- Available EtCO₂ monitoring

Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

Operational Efficiencies

One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
SKU standardization • Inventory reduction • Reduced medical waste
• Reduces complexity for bedside staff



OxyMask™

- Delivers Oxygen flows from 1 to 15+ litres per minute - flush, 24% to 90% FiO₂ (Tyke size 1/4L to 5+, 22% to 65% FiO₂)
- Four sizes available:
 - Adult
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
 - Tyke (6 months - 3 years old, 16 - 33 lbs / 7 - 15 kg)
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



| Code | Description | Mask Size | Tubing Length | Elastic |
|------------|--|-----------|---------------|---------|
| OM-1125-8 | OxyMask™ with 7' tubing | Adult | 7' | Single |
| OM-1125-14 | OxyMask™ with 14' tubing | Adult | 14' | Single |
| OM-1425-8 | OxyMask™ with ear elastics and 7' tubing | Adult | 7' | Double |
| OM-7025-8* | OxyMask™ Non PVC with tubing 7' | Adult | 7' | Single |
| OP-1125-8 | OxyMask™ Plus w 7' tubing | Plus | 7' | Single |
| OP-1125-14 | OxyMask™ Plus w 14' tubing | Plus | 14' | Single |
| OP-1425-8 | OxyMask™ Plus w ear elastics & 7' tubing | Plus | 7' | Double |
| OP-7025-8* | OxyMask™ Plus Non PVC with tubing 7' | Plus | 7' | Single |
| OK-1125-8 | OxyMask™ Kid w 7' tubing | Kid | 7' | Single |
| OK-1125-14 | OxyMask™ Kid w 14' tubing | Kid | 14' | Single |
| OK-7025-8* | OxyMask™ Kid Non PVC with tubing 7' | Kid | 7' | Single |
| OT-1025-8 | OxyMask™ Tyke w 7' tubing | Tyke | 7' | Single |
| OT-7025-8* | OxyMask™ Tyke Non PVC with tubing 7' | Tyke | 7' | Single |

*Non PVC

OxyMask™ EtCO₂

- Oxygen flows from 1 to 12-15 litres per minute
- Delivers 24% to 65% oxygen with End-Tidal CO₂ monitoring
- Three sizes available:
 - Adult
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



| Code | Description | Mask Size | Tubing Length | Elastic |
|---------------|---|-----------|---------------|---------|
| OM-2125-8 | OxyMask™ EtCO ₂ with 7' tubing | Adult | 7' | Single |
| OM-2125-8SLM | OxyMask™ EtCO ₂ , 7' oxygen tubing, 8' sampling line with male connector | Adult | 7' | Single |
| OM-2125-14 | OxyMask™ EtCO ₂ with tubing 14' | Adult | 14' | Single |
| OM-2125-14SLM | OxyMask™ EtCO ₂ , 14' tubing, 15' sampling line with male connector | Adult | 14' | Single |
| OP-2125-8SLM | OxyMask™ EtCO ₂ Plus, 7' oxygen tubing, 8' sampling line with male connector | Plus | 7' | Single |
| OP-2125-8 | OxyMask™ EtCO ₂ Plus with tubing | Plus | 7' | Single |
| OP-2125-14 | OxyMask™ EtCO ₂ Plus dual elastics & with tubing 14' | Plus | 14' | Double |
| OK-2125-8 | OxyMask™ EtCO ₂ Kid with tubing | Kid | 7' | Single |
| OK-2125-8SLM | OxyMask™ EtCO ₂ Kid, 7' oxygen tubing, 8' sampling line with male connector | Kid | 7' | Single |



OxyMulti-Mask™

- Oxygen flows from 1 to 15+ litres per minute - flush
- Delivers 23% to 83% FiO₂ (Kid size 22% to 93% FiO₂)
- Three sizes available:
 - Adult
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
- Delivers oxygen through a standard small bore tube, or using a 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O₂ delivery
- Can be used with, or without high humidity or aerosol therapy
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



| Code | Description | Mask Size | Tubing Length | Elastic |
|--------------|--|-----------|---------------|---------|
| OHH-1025-8 | OxyMulti-Mask™ with adapter & 7' tubing | Adult | 7' | Single |
| OMMK-4025-8 | OxyMulti-Mask™ Adult with 22mm adapter & 7' tubing w/nebulizer | Adult | 7' | Single |
| OPHH-1025-8 | OxyMulti-Mask™ Plus with adapter & 7' tubing | Plus | 7' | Single |
| OPMK-4025-8 | OxyMulti-Mask™ Plus with 22mm adapter & 7' tubing w/nebulizer | Plus | 7' | Single |
| OKHH-1025-8* | OxyUltra Kid™ with adapter & 7' tubing | Kid | 7' | Single |

*OxyUltra Kid™ design not shown

OxyTrach™

- Delivers 23% to 83% FiO₂ on flows from 1 to 15+ litres per minute - flush
- Delivers oxygen through a standard small bore tube using the 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O₂ delivery
- Can be used with, or without high humidity or aerosol therapy
- Comfort strap sold separately
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



| Code | Description | Mask Size | Tubing Length | Elastic |
|-------------|--|-----------|---------------|---------------|
| OAT-1025 | OxyTrach™ 7' tubing | Adult | 7' | Single |
| OAT-1025-CS | OxyTrach™ 7' tubing with comfort strap | Adult | 7' | Comfort Strap |

OxyArm™

- Unique oxygen delivery device with minimal contact and no facial contact
- Ideal for burn units, trauma centers and homecare
- Delivers 24% to 50% oxygen on flows from 1 to 12-15 litres per minute
- Designed for both nose and mouth breathers
- Adapts for left or right sided positioning
- Not made with natural rubber latex
- 25/case (unless otherwise stated)



| Code | Description | Design | Tubing Length |
|----------------|--|------------|---------------|
| OA-DA-1025 | OxyArm™ Disposable Diffuser Arm only | Arm only | none |
| OA-DA-3025 | OxyArm™ Disposable Cannula Arm only | Arm only | none |
| OA-PLUS-1125-8 | OxyArm™ Plus Clear Adjustable Head Band with Diffuser Arm with 7' tubing | Adjustable | 7' |
| OA-PLUS-3125-8 | OxyArm™ Plus Clear Adjustable Head Band with Nasal Cannula Arm & 7" tubing | Adjustable | 7' |
| OA-HB-1125 | OxyArm™ Plus Clear Adjustable Head Band only | Adjustable | none |

