

Sterilization Tapes



- ITSG and JTBSG with STEAM green-to-brown indicator. In natural and blue base material. ITSG also available as latex free version called ITSGL.
- ITE with EO gas purple-to-green indicator.
- ITH with HOT AIR green-to-brown indicator.
- ITP tape strips with VH2O2 blue-to-pink indicator.
- ITP25 tape reel with VH2O2 blue-to-pink indicator.
- Neutral closing tapes with no imprint. In natural and blue base material.

All Indicators are process indicators in accordance with ISO 11140-1 Type 1.

Article	Width	Length	rm/roll	Packing Unit
ITE19	19		50	48
ITH19	19		50	48
ITN19	19		50	48
ITN25	25		50	36
ITN50	50		50	24
ITP20	20	150		200

Article	Width	Length	rm/roll	Packing Unit
ITP25	25		25	24
ITS12G	12		50	72
ITS19G	19		50	48
ITS19GL	19		50	48
ITS25G	25		50	36
ITS25GL	25		50	36
JTB19N	19		50	48
JTB19SG	19		50	48
JTB25N	25		50	36
JTB25SG	25		50	36
TDS25				1

Indicator Strips



Indicator strips are to be placed in each pack of each sterilizer load: Strips changes the colour on correct integration of temperature, time and humidity and concentration of sterilizing agent.

- The AC 250/500 steam strip (type 4) changes colour from white to black
- The AC 6 steam strip (type 6) changes colour from green to brown
- The EO strip (type 5) changes colour from purple to green
- The PL - plasma strip (type 1) changes colour from blue to pink
- The FO strip (type 1) changes colour from dark red to green

in accordance with ISO 11140-1.

Article	Width	Length	Packing Unit
AC250/500	15	200	250
AC6	15	100	250
EO250	15	200	250
FO250/500	15	200	250
PL500	20	150	500

Sterilization Tests



For routine control of autoclave performance:

- Indicator sheets, Daily Test Sheets, for traditional Bowie & Dick type test with a cloth package in accordance with ISO 11140-1 (type 2) and ISO 11140-3 for pre-vacuum steam sterilizers to demonstrate compliance with EN 285. Colour changing from white to black.
- Ready made, compact, single-use Bowie & Dick type test packs: DC20 and DCG2.
- Helix Challenge Test as a product challenge device (PCD) for regular control of hollow loads in autoclaves. Helix consists of the device and 100 indicators strips. Compliance with ISO 11140-1 (type 6) and EN 867-5.

Article	Width	Gusset	Length	Packing Unit
DC20	135	20	105	20
DCG2	100	10	150	20
DT50	210		297	50
HX100	115	23	175	25
HX400	115	23	175	25